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## Continuous Capture Step of a Monoclonal Antibody Downstream Process with Amsphere<sup>™</sup> A3

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There is increased interest in continuous processing for Immunoglobulin G (IgG) purification in bind and elute mode. This is because such continuous chromatography offers higher productivity and cost effectiveness than conventional batch mode. We have developed a continuous chromatography unit operation with a high-capacity monoclonal antibody (mAb) affinity resin, Amsphere<sup>™</sup> A3, and compared its performance with the other resin in the market. Good chromatographic performance was obtained in regard to dynamic binding capacity (DBC) at high flow rate, and impurity clearance such as host cell protein (HCP) and DNA. Also, the results were comparable to those of batch process. These results show that high productivity with good impurity clearance can be expected at continuous capture step using Amsphere A3.

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
Figure 1 Protein A purification process.

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
Figure 2 Breakthrough curve and capacity utilization. Maximum binding capacity of column is described as (a). DBC @10%BT is the area within dotted line at (b). To avoid product loss, sample load amount of batch chromatography is 70% to 90% of DBC. With continuous chromatography, capacity utilization increases described as (c) because second column capture breakthrough product.


Figure 3 Twin-column continuous chromatography.

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Table 1 Productivity simulation conditions

	Continuous process	Batch process
Resin	Amsphere A3	
Column I.D. (cm)	0.5	0.5
Column BH (cm)	5	20
Number of columns	2	1
IgG titer (g/L)	1.5, 4.5	1.5, 4.5
Target Loading amount (g)	15	15
RT (min)	2 min (150 cm/hr) 1 min (300 cm/hr) 0.5 min (600 cm/hr)	4 min (300 cm/hr)
Loading	50%, 70%, 90%	80% of 10% BT

Step	CV*	RT (min)
Wash I	3	CIP: 2 for continuous 8 for batch Other: Same as condition above
Wash II	3	
Wash III	3	
Elution	3	
Strip	3	
EQ1	3	
CIP	3	
EQ2	3	

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Table 2 Purification conditions

	Continuous process	Batch process
Resin	Amsphere A3 50 μm Agarose	Amsphere A3
Column I.D. (cm)	0.5	0.5
Column BH (cm)	5	10
Number of columns	2	1
RT (min)	0.75 min (400 cm/hr) 0.5 min (600 cm/hr)	1.5 min (400 cm/hr)
Loading	70 %BT	80 % of 10 % BT
Number of cycles	4	1

Step	Buffer	CV
Load	Bevacizumab BS, Titer 1.6 g/L	-
Wash I	20 mM Sodium Phosphate, pH7.5	3
Wash II	20 mM Sodium Phosphate, 1 M NaCl, pH7.5	3
Wash III	20 mM Sodium Phosphate, pH7.5	3
Elution	100 mM Sodium Acetate, pH3.3	3
Strip	100 mM Acetic Acid, pH2.8	3
EQ1	20 mM Sodium Phosphate, pH7.5	3
CIP	0.5 M NaOH	3
EQ2	20 mM Sodium Phosphate, pH7.5	3

\*CIP: RT (Residence Time) = 2 min, 150 cm/hr (continuous process), RT = 4 min, 150 cm/hr (batch process), Other: RT is same as condition above.


Figure 4 DBC comparison between Amsphere A3 and Agarose resin.


Figure 5 Productivity of Amsphere A3 with various conditions.

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Table 3 Recovery and measured productivity

Resin	Continuous process			Batch process
	Amsphere A3		50 μm Agarose	Amsphere A3
RT (min)	0.75 min (400 cm/hr)	0.5 min (600 cm/hr)	0.75 min (400 cm/hr)	1.5 min (400 cm/hr)
Loading	70 % BT	70 % BT	70 % BT	80 % of 10 % BT
Binding Capacity (g/L-resin)	65	49	71	34
Calculated Productivity (g/(L·hr))	30	35	29	27
Recovery (%)	98	99	95	97
Measured Productivity (g/(L·hr))	30	35	28	26


Figure 6 Impurity clearance.

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BioProcess International Conference & Exhibition 2017, 2017/9/25-28

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